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


PROFILE DESIGNER
i physio®


NOW AVAILABLE ON MULTI UNIT ABUTMENTS !


2 EVOLUTIONS FOR YOUR MULTIPLE IMPLANTS CASES

- New design inspired from iPhysio for our « Multi Unit » type abutments
- New range of iPhysio for « Multi Unit » type abutments

DISCOVER THE BENEFITS OF IPHYSIO FOR YOUR TETRA AND PLURAL ABUTMENT RESTORATIONS

- **iPhysio emergence profile** seamlessly continuous with the TETRA abutment
- **Textured surface** with a **zirconium nitride coating**
- **One component, 3 functions :**
 - Protects the TETRA abutment and promotes anatomical soft tissue healing
 - Functions as a **scan body** for digital impressions
 - Provides support for provisional restorations

THE GOLD STANDARD IN IMPLANTOLOGY

- **No sharp edges**, providing greater patient comfort.
- **Simplified continuous scanning** from gingiva to iPhysio for a smooth and highly accurate digital impression.
- **Scanning time reduced by up to threefold**—no more long and tedious scanning of scan bodies protruding through the gingiva. Only the **occlusal surface of the iPhysio** needs to be scanned; there is no need to capture the transmucosal emergence profile beneath the iPhysio.
- **Outstanding scanning accuracy** for multiple-unit cases thanks to **iPhysio technology** (micro-textured zirconium nitride-coated surface) and the excellent scanability of the iPhysio geometry.
- **Clear identification of individual iPhysio heads**, eliminating the risk of confusion with adjacent scan bodies, thanks to their low-profile position close to the gingival margin.

INDICATIONS

- **Small bridges from 2 to 6 elements**
- **on 2 to 4 implants**
- Do not confuse with full arch : use our new iFlag Reverso (available early 2027).
- For single restorations : a new version of our tetra abutment is coming soon.

PROTOCOL

1 . Implants placement

- Avoid placing implants in perfect parallel alignment to prevent the "three-legged stool with aligned legs" effect, which compromises stability.
- Prioritize implant placement where sufficient bone is available, and use **angled TETRA abutments** when necessary to correct prosthetic angulation.
- Perfect implant alignment is an outdated requirement that was primarily driven by the limitations of conventional analog impression techniques.
- With digital workflows, it is preferable to prioritize optimal load distribution and overall prosthetic stability.

2 . Placing the Tetra abutment

- The height of the **TETRA abutment** depends on the depth of implant placement and the gingival height.
- It should be positioned **1 mm subgingivally** to respect the iPhysio concept, as its convex profile is designed to support the soft tissues (**0.5 mm** when 1 mm is not achievable).
- For **maxillary central incisors**, a **narrower PLURAL abutment** should be selected.
- **Tightening torque: 35 N·cm**

3 . Selecting and placing the iPhysio

- In most clinical situations, the recommended choice is **commercial height 2**.
- Please note that this designation corresponds to the height of the former protective cap (**actual overall height: 4.2 mm**). It is suitable for almost all cases except when the implant is placed deeply subcrestally (approximately **2 mm**) in combination with thick soft tissue (**6 mm**).
- The **actual height** of the **commercial height 4 iPhysio** is **6 mm**.
- **Tightening torque: 10 N·cm**
- To stabilize the rotational iPhysio during tightening, use **Führer forceps**.

4. Healing

- Allow a healing period of **at least 3 months**, depending on the clinical site.

5. Digital impression on iPhysio

- The accuracy of the digital impression has been validated for the indicated clinical applications.

6. Provisional restoration

- The prosthetic axes should be sufficiently parallel to facilitate removal of the provisional restoration.
- When implant angulations need to be corrected, use **angled TETRA abutments**.
- The retention of the provisional restoration is intentionally high (**factory-calibrated to resist approximately 14 N of pull-off force**) to minimize the risk of bridge dislodgement.
- If lower retention is desired, reduce it by milling **2 of the 6 retention fins** that provide the snap-fit mechanism.

COMPATIBLE WITH ALL MULTI-UNIT TYPE ABUTMENTS

PROFILE DESIGNER
i physio[®]